



**EASTERN REGIONAL RESEARCH CENTER  
AGRICULTURAL RESEARCH SERVICE  
UNITED STATES DEPARTMENT OF AGRICULTURE  
600 E. MERMAID LANE  
WYNDMOOR, PA 19038  
(215) 233-6400**

**Title:** *Cyclospora cayetanensis*: Research Methodology and Control by Irradiation

**Author(s):** B.A. Niemira, and D.W. Thayer

**Citation:** In "Use of Irradiation to Ensure Hygienic Quality of Fresh, Pre-Cut and Vegetables and Other Minimally Processed Food of Plant Origin" Report of First FAO/IAEA Meeting 5-9, November 2001, Rio De Janeiro, 7 pp. (2001) 1-7

**Number:** 7439

**Please Note:**

This article was written and prepared by U.S. Government employees on official time, and is therefore in the public domain.

Our on-line publications are scanned and captured using Adobe Acrobat. During the capture process some errors may occur. Please contact William Damert, [wdamert@arserrc.gov](mailto:wdamert@arserrc.gov) if you notice any errors in this publication.

1 

# ***Cyclospora cayetanensis*:** **Research methodology and control by** **irradiation**

**Brendan A. Niemira and**  
**Donald W. Thayer**

Food Safety Intervention Technologies Research Unit  
U.S. Dept. Agriculture, Agricultural Research Service  
Eastern Regional Research Center  
Wyndmoor, Pennsylvania

2 

## ***Cyclospora cayetanensis*:** **introduction and epidemiology**

- Coccidian parasite; spherical oocysts ~9 nm in diameter
- Reports of recovery of *C. cayetanensis* from Old World primates, but humans may be the only true hosts
- Disease progression:
  - ingestion of a sporulated oocyst containing 4 sporozoites
  - sporozoites penetrate epithelial lining of intestine
  - increasing damage to host cells concurrent with parasite's asexual division, mobility and sexual reproduction
  - new oocysts sloughed into feces, along with damaged epithelium

3 

4 

## ***C. cayetanensis*:** **introduction and epidemiology**

- Clinical symptoms
  - signs start 1-7 days post-ingestion of oocysts
  - mild infections may produce little or no outward signs
  - more extensive infections may show prolonged watery diarrhea, abdominal cramping, weight loss, anorexia, vomiting and/or fever
- Extensive treatment with sulfa-type drugs can be effective

5 

## ***C. cayetanensis*:** **introduction and epidemiology**

- Infection historically associated with poor sanitation
  - pediatric gastroenteritis in developing nations

- adult infection in travelers from industrialized nations visiting loci of infection
- 1996: 1465 cases in USA, Canada associated with Guatemalan raspberries (Herwaldt and Ackers, N Engl J Med 1997 May 29;336(22):1548-56)
- April-June 1997: 1012 cases in USA, Canada, Guatemalan raspberries again implicated (Herwaldt and Beach, Ann Intern Med 1999 Feb 2;130(3):210-20)
- May 1997: Guatemala voluntarily suspends export of fresh raspberries, outbreak ends

## 6 *C. cayetanensis*:

### introduction and epidemiology

- Outbreaks in North America occur more frequently in spring, early summer
- Relation to time of peak importation of fresh fruits and vegetables suggested
- Long distance transport of commodities, multiple handling stations increase possibility of contamination between farm and table
- Washing foods before consumption is a key step, but it does not remove 100% of oocysts (Ortega, Y.R. et al. 1997 Am. J. Trop. Med. Hyg. 57: 683-686)

## 7 *C. cayetanensis*: detection methodology

- US-FDA BAM: Concentration and Preparation of Cyclospora from Berries for the Polymerase Chain Reaction (PCR) and Microscopy ([vm.cfsan.fda.gov/~ebam/bam-19a.html](http://vm.cfsan.fda.gov/~ebam/bam-19a.html))
- Epi-illuminated Fluorescence Microscope; UV 1A filter block (Excitation Filter, EX 365/10; Dichroic mirror, DM 400; Barrier Filter, BA-400; or equivalent)
- PCR: buffers, reagents, thermocycler, gels, film (protocols should be validated for each material/isolate)

## 8

- Wash method for fresh produce (berries, lettuce, etc.) or puree
- Stomacher bag, 250 ml of deionized water
- 250-500 g of produce
  - berries +/- 1 berry; USE ONLY INTACT BERRIES
  - juice from cut or broken berries may be inhibitory to PCR
  - debris may interfere with microscopy
- 250 g of puree (juice and debris are unavoidable)

## 9

- Invert and agitate gently. Avoid damaging produce. Centrifuge 1500 X g for 10 minutes
- Decant supernatant, retaining 1 ml of supernatant and to resuspend pellet fraction.
- Appropriate aliquots of resuspended pellet will be removed from this tube for microscopy and PCR.
- Store the remaining portion at 4°C for up to one month. After one month, dilute the remaining material with an equal volume of 2.5% potassium dichromate, mix, and store at 4°C.

## 10

- Prepare wet mount of 10 µl sediment.
- View under UV light at 400 X. Cyclospora oocysts fluoresce cobalt blue with UV-1A filter, blue-green with broader UV illumination.
- Confirm cyst size of 8-10 µm with microscope reticle, compare presumptive oocysts to those in a known standard.
- Switch from epi-fluorescence microscopy to bright field or DIC to confirm internal structures
- Seal slides, document positive samples with photographs taken at multiple planes.

11 

- 100 µl of produce sediment. Thaw, microcentrifuge. 14,000 RPM (15,800 X g) for 3 min and discard supernatant.
- Wash pellet once with 100 µl TE buffer, centrifuge at 14,000 RPM for 3 min. Discard supernatant and resuspend in 100 µl TE
- Complete 3 freeze/thaw cycles, each 2 min in liquid nitrogen or a dry ice-ethanol bath followed by 2 min in a 98°C water bath. Add 0.1 ± 0.02 g glass beads to extract.

12 

- Vortex for 5 min, chill on ice for 5 min.
- Centrifuge sample extract at 14,000 RPM (15,800 X g) for 3 min. Transfer supernatant to new sterile microcentrifuge tube.
- This extract can be stored frozen (-20°C) until ready for the PCR analysis, or if needed as reserve in case of template inhibition problems.
- Combine 20 µl sample extract and 2 µl freshly prepared Non-fat Milk Solution in 1 ml sterile deionized water. The entire 22 µl will be used as template in a 100 µl PCR amplification.

13 

- PCR primers: CYCF1E, CYCR2B, CYCF3E and CYCR4B (Relman et. al., J. Infect. Dis. 173:440-445, 1996)

14 

- Two sequential PCR runs
- #1 uses CYCF1E and CYCR2B primers, probing the template DNA extracted from the sample
- #2 uses CYCF3E and CYCR4B primers, probing the template DNA from Step #1
- Agarose/ethidium bromide gel electrophoresis, UV transillumination. Predicted size of DNA marker after F1E/R2B, is 651 bp; after F3E/R4B is 308 bp.
- ★ Amplified product after the first round may not be visible; only product from the second round of PCR should be electrophoresed.

15 

- PCR product of 308 bp after the second PCR round a presumptive positive for Cyclospora or Eimeria.

- Digestion of the amplification product with MnlI distinguishes *Cyclospora* from *Eimeria*
- Prepare separate restriction digests
  - each presumptive positive PCR amplification product
  - amplification products from control *Cyclospora cayetanensis* and *Eimeria tenella* strains.
- Agarose/ethidium bromide gel electrophoresis

16 ☐

- Microscopy is a fundamental tool for detection
- PCR will reliably indicate the presence of *C. cayetanensis*, and can discriminate between *Cyclospora* and *Eimeria*
- For both methods, comparison with known standards is essential
- ★ PCR will amplify DNA from dead cells - not a diagnostic tool to evaluate post-irradiation survival or infectivity. Only indicates presence/absence.

#### 17 ☐ *C. cayetanensis*: ionizing radiation

- Multiple studies have been completed on *C. cayetanensis* (Dubey, Thayer, Speer, and Shen, 1998)
- Typically trials were conducted with 50,000 sporulated or 50,000 unsporulated oocysts
- Excystation occurred at doses through 0.5 kGy
- Sporulation occurred at doses up to 1.0 kGy though there was a dose dependent delay, morphological changes, and reduction in sporulation

#### 18 ☐ *C. cayetanensis*: treatment validation

- No animal model exists; difficult to validate inactivation
- Molecular techniques may provide some additional tools, but validation of inactivation treatments requires human experimentation
- Severity of symptoms, frequency of mortality hinders efforts to recruit volunteers
- A good model pathogen would provide results that would reliably predict the behavior of irradiated *C. cayetanensis*

#### 19 ☐ *Toxoplasma gondii* oocyst as model for coccidian parasite irradiation

- Low host specificity, infects many mammals, including human. Good animal model (mouse) readily available.
- Infection generally asymptomatic in adults; low risk of hepatitis, blindness, spontaneous abortion
- Transmitted by fecal, oral and meat
- Very efficient bioassay — 1 oocyst can be fatal in test subject.
- Mortality in mice correlated with dose exposure

#### 20 ☐ *T. gondii*: methodology

- Dubey, Jenkins, Thayer 1996. J. Parasitol. 82:724-727
- Oocysts obtained by sugar flotation of feces from cats fed tissues cysts of VEG strain.
- Oocysts sporulated in 2% H<sub>2</sub>SO<sub>4</sub> at room temperature, 7 days. Stored in H<sub>2</sub>SO<sub>4</sub> at 4C until use.
- ~10<sup>6</sup> sporulated oocysts neutralized in 3% NaOH, centrifuged, resuspended in 1 ml 0.9% NaCl
- Gamma irradiation was performed at room temperature
- Treated oocysts were inoculated orally into mice

#### 21 ☐ *T. gondii*: methodology

- Mice that died after inoculation were examined for tissue cysts
- Surviving mice were bled periodically to isolate antibodies to *T. gondii*
- All mice were sacrificed after 2 months. A portion of the mouse tissue was fed to 2 *T. gondii*-free cats.
- The cats were bled to isolate antibodies, and sacrificed after 25 days to identify tissue cysts

#### 22 ☐ *T. gondii* radiation biology

- Most mice fed non-irradiated oocysts or oocysts irradiated to 0.1 kGy or less died of acute toxoplasmosis
- At doses 0.2 and above, no tissue cysts were found.
- Antibodies to *T. gondii* were found in 11 of 40 mice fed oocysts irradiated to 0.2-0.5 kGy
- Cats fed tissue from seropositive mice did not shed *T. gondii* oocysts in feces, did not develop antibodies to *T. gondii* and did not develop tissue cysts

#### 23 ☐ *T. gondii*: strain radiation sensitivity

- Dubey and Thayer, J. Parasitol. 80:764-767 (1994)
- Mice were inoculated with 95 new strains of *T. gondii*, isolated from pigs. Also inoculated with 10 existing laboratory strains
- Brain tissue vacuum packed, irradiated to 0.0 (control), 0.1, 0.2, 0.3, 0.4, 0.5, 0.6, 0.7, and 0.9 kGy
- Treated samples were bioassayed in mice, cats or both
- Tissue cysts of all strains were rendered nonviable at 0.4 kGy
- 0.25 kGy had identical efficacy when delivered at -4, 0, 4, 8, 12 or 16C. No effect of irradiation temperature.

#### 24 ☐ *T. gondii*: oocysts irradiated before sporulation

- Unsporulated oocysts irradiated
  - Exp't 1: 0.0 (control), 0.2, 0.4, 0.6 and 0.8 kGy
  - Exp't 2: 0.0 (control), 3.5 and 4.0 kGy
- Incubated in 2% H<sub>2</sub>SO<sub>4</sub> at room temperature for 30 (#1) or 12 (#2) days to allow sporulation
- Microscopy to confirm germination, fed to mice

- Development of cysts evaluated by bioassay
- 25 ☐ *T. gondii*: oocysts irradiated before sporulation
- \* Antibodies to *T. gondii* not found in these mice. No evidence of disease demonstrated by bioassay
  - (Dubey et al. 1998 Int. J Parasitology 28:369-375)
- 26 ☐ *T. gondii*: effect of sporulation stage
- \* Oocysts were irradiated before sporulation (2 sporocysts formed) and then bioassayed in mice.
  - \*\* No. of mice *T. gondii* positive of no. of mice inoculated.
- 27 ☐
- Oocysts sporulated before irradiation
  - Raspberries either sprayed via atomiser ( $10^6$  oocysts) or injected with suspension ( $10^4$  oocysts)
  - Both sets irradiated to 0.4 kGy
  - Berries homogenized in saline, pelletized. Sediment filtered (200 and 90  $\mu$ m screen) and repelletized.
  - Sediment fed to mice. Bled for antibody isolation and sacrificial bioassay.
- 28 ☐ *T. gondii*: oocysts irradiated on raspberries
- 29 ☐ *T. gondii*: conclusions
- *T. gondii* sporulates after doses up to 1.0 kGy, but is nonviable and unable to cause disease after doses of 0.4 kGy
  - *T. gondii* oocysts irradiated at  $\geq 0.4$  kGy can excyst but do not multiply in the host
  - Sporozoites from irradiated oocysts can induce antibody formation in the host, even in the absence of histopathological signs
  - *In vitro* measurement of sporulation alone insufficient to properly evaluate efficacy of irradiation
- 30 ☐ *T. gondii*: conclusions
- Biology and radiation sensitivity make *T. gondii* a good model for *C. cayetanensis*
  - Requirement for animal model to validate efficacy of irradiation, as *in vitro* observations do not reliably indicate ability to cause disease
  - Elimination/inactivation of *T. gondii* at relatively low doses,  $\sim 0.4$  kGy
  - This dose level is tolerable and beneficial for many fruits and vegetables, including raspberries
  - US-FDA regulations: max. of 1.0 kGy for disinfestation, delay of maturation of produce. Proposal to increase limit is under review.
- 31 ☐ Research questions
- Effect of different fruit and/or vegetable substrates on radiation sensitivity of *C. cayetanensis* or *T. gondii*

- Effect of isolate and/or strain variation, in combination with changes in substrate, processing conditions (temperature, MAP, etc.)
- Radiation processing of other coccidians (e.g. *Cryptosporidium*)
- Further validation and acceptance of *T. gondii* as a model for *C. cayetanensis* by regulatory bodies, and interaction of different regulatory agencies in exporting/importing countries

## 32 Electronic reference resources

- US-FDA: Concentration and Preparation of Cyclospora from Berries for the Polymerase Chain Reaction (pcr) and Microscopy
- US-CDC: *C. cayetanensis* FAQ, epidemiology links  
([www.cdc.gov/ncidod/dpd/parasites/cyclospora/factsht\\_cyclospora.htm](http://www.cdc.gov/ncidod/dpd/parasites/cyclospora/factsht_cyclospora.htm))
- Kansas State University Parasitology Laboratory  
([www.ksu.edu/parasitology/cyclospora/cyclospora.html](http://www.ksu.edu/parasitology/cyclospora/cyclospora.html))
- USDA: Food Irradiation (overview, including discussion of *C. cayetanensis*)  
([www.ars.usda.gov/is/pr/1997/971210.htm](http://www.ars.usda.gov/is/pr/1997/971210.htm))